

Technical Data Sheets



Product Name: COUCH COVER ROLL LAMINATED

Premium Standard Laminated Paper + PE

50 CM X 50 M

Product Code: M951786R

Paper Specifications:				
Specifications	Standard	Unit	Data 1	Data 2
Composition	100% virgin pulp	-	-	
Brightness	UNI EN ISO 12625-3	%	86	+/- 5%
Grammage	UNI EN ISO 12625-6	gr/m ²	18Gr/m2	+/- 7%
Ply		Ply		
MD Tensile strength	UNI EN ISO 12625-4	grf/25mm	>1650	+/- 10%
CD Tensile strength	UNI EN ISO 12625-4	grf/25mm	>1100	+/- 10%
MD Wet Tensile strength	UNI EN ISO 12625-5	grf/25mm	>300	+/- 5%
CD Wet Tensile strength	UNI EN ISO 12625-5	grf/25mm	-	+/- 5%
MD Elongation	UNI EN ISO 12625-4	%	20-26	+/- 5%

PE Specifications:				
Specifications	Standard	Unit	Data 1	Data 2
Composition:	100% High Density	-	-	
Colour		-	White	
Grammage	UNI EN ISO 12625-6	g/m ²	8	+/- 5%

Roll Specifications:			
Specifications	Unit	Data 1	Data 2
Roll width:	cm	50	+/- 1%
Roll length:	m	50	+/- 2%
Roll weight:	g	726	+/- 5%
Roll diameter:	cm	11	+/- 5%
Core diameter:	cm	4,3	+/- 1%
Sheet length:	cm	50	+/- 1%
Sheet number:	No	100	+/- 1%

Packaging Specifications:				
Sütun1	Dimension	Quantity	Weight	Data 2
Roll:	50 CM X 50 M	1 Roll	726 Grams	+/- 5%
Carton:	400x535x290mm	12 Rolls	8.71 Kg	+/- 5%
Pallet:	120 x 80 x 224	28 Cartons	243.88 Kg	+/- 5%

ADDRESS
Ergazi Mahallesi 1819 Cadde No:3D
(TEKNOSİTE) Yenimahalle / ANKARA /
TURKIYE

Tel: +90 312 255 60 03 (pbx)
Fax: +90 312 255 60 23

